

Connecting your hearing aids to the ReSound 3D app for remote control access:

1. Begin by downloading the Resound Smart 3D app and tapping on it. You can access this from your app store (iphone), or play store (android).
2. You'll see the "welcome screen", and you need to then press "get started".
3. Next, some terms and conditions will appear, which you should read and then tap "agree" if you're happy to continue.
4. Tap "OK" when it asks if the app can use your phone's Bluetooth function.
5. Then open the battery doors of your hearing aids and close them.
6. You might have to repeat this step if indicated on the app.
7. Once the app says "connected", press "continue".
8. You can then choose to "allow" or "don't allow" notifications – a personal preference.
9. Press "use location" when prompted.
10. Your should now be ready to use your ReSound Smart 3D app.